

Ex vivo catheter-based imaging of coronary atherosclerosis using multimodality OCT and NIRAF excited at 633 nm

March 19, 2015

While optical coherence tomography (OCT) has been shown to be capable of imaging coronary plaque microstructure, additional chemical/molecular information may be needed in order to determine which lesions are at risk of causing an acute coronary event. In this study, we used a recently developed imaging system and double-clad fiber (DCF) catheter capable of simultaneously acquiring both OCT and red excited near-infrared autofluorescence (NIRAF) images (excitation: 633 nm, emission: 680nm to 900nm). We found that NIRAF is elevated in lesions that contain necrotic core – a feature that is critical for vulnerable plaque diagnosis and that is not readily discriminated by OCT alone.

Authors : Hao Wang, Joseph A. Gardecki, Giovanni J. Ughi, Paulino Vacas Jacques, Ehsan Hamidi, Guillermo J. Tearney

Please visit Biomedical Optics Express for access to the full paper

<https://www.osapublishing.org/boe/abstract.cfm?uri=boe-6-4-1363>